



HP Confidential. For use by HP or Partner with Customers under HP CDA only.

Please refer to the HP Healthcare Edition HC270cr Clinical Review display user guide for instructions on cleaning this product. Here are additional instructions for cleaning with specific germicidal wipes.

- First turn off the computer and display to prevent electric shock or damage to components.
 - Disconnect AC power
 - Disconnect all powered external devices
- If there is any visible soiling, wipe the exterior of the display with a soft, water-dampened cloth to remove the soil as needed.
 - The cloth should be of dry microfiber or a chamois (static-free cloth without oil), or static-free cloth wipes.
 - The cloth should be moist but not wet. Water dripping into the ventilation or other points of ingress can cause damage.
 - Please do not use fibrous materials, such as paper towels, as these can scratch the display. Over time, dirt and cleaning agents can get trapped in the scratches.
- Next, use any of the following germicidal wipes* to clean all exposed surfaces of your display. Remember to start with the display and finish with any flexible cables. (Refer to the directions for use provided by the manufacturer of the wipes).
 - **Sani-Cloth® HB**
 - Benzyl-C12-18-alkyldimethyl ammonium chlorides: <0.1%
 - Quaternary ammonium compounds, C12-14-alkyl[(ethylphenyl)methyl]dimethyl, chlorides: <0.1%
 - **Sani-Cloth® Plus**
 - Isopropanol: 10-20%
 - 2-Butoxyethanol: <5%
 - Quaternary ammonium compounds, C12-18-alkyl [(ethylphenyl) methyl] dimethyl, chlorides: <0.5%
 - Benzyl-C12-18-alkyldimethyl ammonium chlorides: <0.5%
 - **Sani-Cloth® AF3**
 - Quaternary ammonium compounds, C12-18-alkyl [(ethylphenyl) methyl] dimethyl, chlorides: <0.5%
 - Benzyl-C12-18-alkyldimethyl ammonium chlorides: <0.5%
 - **Super Sani-Cloth® Germicidal Disposable Wipes**
 - Isopropyl alcohol: 55.0%
 - Alkyl dimethyl benzyl ammonium chlorides: 0.250%
 - Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride: 0.250%
 - **CaviWipes®**
 - Isopropanol: 10-20%
 - Ethylene Glycol Monobutyl Ether (2-Butoxyethanol): 1-5%
 - Diisobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride: 0.1-0.5%
 - **Clorox Healthcare® Bleach Germicidal Wipes**
 - Sodium hypochlorite 0.1%–1.0%
 - **Oxivir® Tb Wipes**
 - Hydrogen Peroxide: 0.1-1.5%
 - **Windex® Electronic Wipes**
 - Cellulose: 10-30%
 - Ethyl alcohol: 0.10-1.00%
 - **PDI Easy Screen® Cleaning Wipe**
 - Isopropanol: 60-70%
 - Water: 30-40%
- Allow the unit to air-dry before use

Please do not use cleaners that contain any petroleum-based materials such as benzene or thinner. These may damage the display.

Sincerely,

© COPYRIGHT 2018 * DISCLAIMER: Solutions tested as of August 18, 2017; chemical composition is subject to change. Repeated use of these germicidal wipes may cause some cosmetic changes to the product. 4AA7-3877ENW

hp.com/go/healthcare